



REACTHEALTH

# BILEVEL COMPARISON CHART

**MANUFACTURER**

*React Health*



*Philips*



*ResMed*



BRAND	G3 BPAP 25A	DreamStation® BiPAP Pro	AirCurve® 10 VAuto
<b>Dimensions</b>	5.70 x 10.48 x 4.49 in.	11.7 x 7.6 x 3.3 in.	4.57 X 10.04 X 5.91 in.
<b>Weight</b>	3.74 lbs.	4.37 lbs.	3.52 lbs.
<b>Humidification</b>	Standard	Fixed, Adaptive, Heated Tube modes	Standard, Integrated, Fixed or Adaptive
<b>Heated Tubing</b>	Standard	Optional	Optional
<b>Connectivity</b>	iCode connect, QR Code, SD Card, Integrated Cellular Modem	Integrated Bluetooth, optional Wi-Fi & cellular modems	Wi-Fi, SD Card
<b>SD Data Capacity</b>	Minimum 12 months	> 1 year on-board: <6 months	N/A
<b>Filters</b>	Reusable inlet filter	Reusable pollen filter, disposable ultra-fine filter	Disposable Antibacterial Filter
<b>Modes</b>	CPAP, AutoCPAP, S, AutoS	Fixed CPAP, Fixed BiPAP	CPAP, BiLevel, Spontaneous
<b>Rise Time</b>	150-400 ms	100-400 ms	150-900 ms
<b>Ramp Time</b>	0 to 60 minutes	5 to 45 minutes (5 minute increments)	0 to 45 minutes, user adjustable
<b>Ventilation Pressure</b>	4 to 25 cm H2O	4 to 25 cm H2O	3-25cm H2O
<b>Weight Restrictions</b>	Over 66 lbs	Over 66 lbs.	Over 66 lbs.

Dreamstation® is a registered trademark of Koninklijke Philips N.V. (Royal Philips)  
AirCurve® is a registered trademark of ResMed family of companies.

\*React Health is formerly known as 3B Medical, Inc.