

PAP COMPARISON CHART

MANUFACTURER

React Health

React Health

ResMed

Philips

Resvent



| BRAND | LUNA II | LUNA G3 | AIRSENSE® 11 | DREAMSTATION® 2 | iBreeze |
|-----------------------------|-----------------------------------|-----------------------------------|---------------------------------------|--------------------------|--------------|
| MODES | CPAP, Auto CPAP | CPAP, Auto CPAP | AutoSet, CPAP, Elite | Auto | Auto CPAP |
| Pressure Range | 4-20 | 4-20 | 4-20 | 4-20 | 4-20 |
| Exhalation Relief | Yes (RESlex) | Yes (RESlex) | Yes (EPR) | Yes (Flex) | Yes (IPR) |
| Remote Setting Changes | Yes | Yes | Yes | Yes | No |
| Integrated Humidifier | Yes | Yes | Yes | Yes | Yes |
| Sound Level dBa | 28dBa | 26dBa | 27dBa | 28dBa | 28dBa NA |
| Compliance Reporting | Yes | Yes | Yes | Yes | Only SD Card |
| Phone App | Yes | Yes | Yes | Yes | Yes |
| Connectivity | Cellular, Wi-Fi, QR Code, SD Card | Cellular, Wi-Fi, QR Code, SD Card | Cellular, SD Card | Cellular, Wi-Fi, SD Card | NA |
| Heated Tubing | Optional | Standard | Optional (subject to availability) | Optional | Optional |
| On Screen Resupply Reminder | No | Yes | No | No | No |