

# Introducing the Luna Family of PAP Devices from React Health



Luna G3 CPAP (LG3500) and APAP (LG3600)



Luna G3 Bilevel 25A (LG3700) and Bilevel ST 30Vt (LG3800)

## LUNA FEATURES:

### The Luna Family of PAP Devices offer:

- Reslex™ pressure relief on exhalation feature for increased patient comfort and compliance
- Leak compensation to maintain patient's pressure and breathing synchronization
- Simple release water chamber for easier removal and cleaning
- Large screens and easy to navigate menus for clinicians and patients
- Integrated heated humidification
- Event detection for apneas, hypopneas, snoring and central apneas
- Bilevel ST model offers Target Volume
- Heated Tubing standard on all Luna G3 Models
- **FDA Approved**

### Luna PAP Devices Make Capturing Compliance Easy,

#### 5 separate ways to retrieve patient data:



Cellular



Wi-Fi



SD Card



QR Code Scanning



iCode Number

Patient data uploads to *iCodeConnect*, React Health's Cloud based platform for viewing patient data.



REACTHEALTH

# PAP COMPARISON CHART

## MANUFACTURER

React Health

React Health

ResMed

Philips

Resvent



BRAND	LUNA II	LUNA G3	AIRSENSE® 11	DREAMSTATION® 2	iBreeze
MODES	CPAP, Auto CPAP	CPAP, Auto CPAP	AutoSet, CPAP, Elite	Auto	Auto CPAP
Pressure Range	4-20	4-20	4-20	4-20	4-20
Exhalation Relief	Yes (RESlex)	Yes (RESlex)	Yes (EPR)	Yes (Flex)	Yes (IPR)
Remote Setting Changes	Yes	Yes	Yes	Yes	No
Integrated Humidifier	Yes	Yes	Yes	Yes	Yes
Sound Level dBa	28dBa	26dBa	27dBa	28dBa	28dBa NA
Compliance Reporting	Yes	Yes	Yes	Yes	Only SD Card
Phone App	Yes	Yes	Yes	Yes	Yes
Connectivity	Cellular, Wi-Fi, QR Code, SD Card	Cellular, Wi-Fi, QR Code, SD Card	Cellular, SD Card	Cellular, Wi-Fi, SD Card	NA
Heated Tubing	Optional	Standard	Optional (subject to availability)	Optional	Optional
On Screen Resupply Reminder	No	Yes	No	No	No