



REACTHEALTH

BILEVEL COMPARISON CHART

MANUFACTURER

React Health



Philips



ResMed



BRAND	G3 BPAP 25A	DreamStation® BiPAP Pro	AirCurve® 10 VAuto
Dimensions	5.70 x 10.48 x 4.49 in.	11.7 x 7.6 x 3.3 in.	4.57 X 10.04 X 5.91 in.
Weight	3.74 lbs.	4.37 lbs.	3.52 lbs.
Humidification	Standard	Fixed, Adaptive, Heated Tube modes	Standard, Integrated, Fixed or Adaptive
Heated Tubing	Standard	Optional	Optional
Connectivity	iCode connect, QR Code, SD Card, Integrated Cellular Modem	Integrated Bluetooth, optional Wi-Fi & cellular modems	Wi-Fi, SD Card
SD Data Capacity	Minimum 12 months	> 1 year on-board: <6 months	N/A
Filters	Reusable inlet filter	Reusable pollen filter, disposable ultra-fine filter	Disposable Antibacterial Filter
Modes	CPAP, AutoCPAP, S, AutoS	Fixed CPAP, Fixed BiPAP	CPAP, BiLevel, Spontaneous
Rise Time	150-400 ms	100-400 ms	150-900 ms
Ramp Time	0 to 60 minutes	5 to 45 minutes (5 minute increments)	0 to 45 minutes, user adjustable
Ventilation Pressure	4 to 25 cm H2O	4 to 25 cm H2O	3-25cm H2O
Weight Restrictions	Over 66 lbs	Over 66 lbs.	Over 66 lbs.

Dreamstation® is a registered trademark of Koninklijke Philips N.V. (Royal Philips)
 AirCurve® is a registered trademark of ResMed family of companies.

*React Health is formerly known as 3B Medical, Inc.